PUSH MED ANKLE BRACE

AEQUI FLEX

INTRODUCTION

Push med products provide effective solutions for the treatment and prevention of injuries of the locomotor apparatus. Approved by physicians and scientists.

Your new ankle brace is a Push med product, characterised by a high level of effectiveness, comfort and quality. The Push med Aequi flex increases ankle stability and helps to control pain. This increases physical strength. Mechanical support, adjustable pressure and good position retention constitute the basis of its performance.

Push med products are innovative and are manufactured for you with great care. This ensures maximum security and more freedom of movement. To promote satisfied use.

Read the product information carefully before you use the Push Aequi flex. If you have any questions, please consult your reseller or send an e-mail to info@push.eu

PRODUCT CHARACTERISTICS

The Push Aequi flex is an effective ankle brace, which strikes a balance between efficient mechanical support of an unstable ankle and good wear comfort. The rigid support is situated exclusively on the medial (inner) side of the ankle and has a cavity for the malleolus. The lateral side is made from comfortable and elastic foam, which easily adapts to the individual anatomy. The bracing material is constructed in such a way that it is rigid in the centre, becoming softer towards the edges to ensure a good fit and wearer comfort. The brace is fastened by means of a diagonal strap. The elastic straps apply pressure around the ankle. The Push Aequi flex can be used in the functional treatment of ankle injuries and as a long term support in cases of chronic ankle instability. The inversion/eversion motion (inward / outward movement of the foot) is restricted while plantar/dorsal flexion (bending and stretching) is allowed, for adequate freedom. This means that no inhibition is felt during normal walking and running. The numbers on the straps ensure easy fitting. The Push Aequi flex is worn over the sock, and optimal stabilisation can only be achieved if it is worn with a closed (lace-up) shoe.

INDICATIONS

- Treatment of acute lateral ankle ligament lesions
- Nachbehandlungen von konservativen oder operierten Knöcheltraumata
- Follow-up treatment of conservative or operated ankle fractures.
- Serious (residual) instability
- Arthrosis with instability
- Secondary prevention of ankle ligament lesions

There are no known disorders where the brace cannot or should not be worn. In the case of poor blood circulation in the limbs, a brace may cause hyperaemia. Avoid applying the brace too lightly in the beginning.

Pain is always a sign of warning. If the pain persists or the complaint deteriorates, it is advisable to consult your doctor or therapist.

FITTING

The Velcro fastenings are marked with symbols. These are to help you fasten the Velcro in the correct sequence. First  on , then  on  etc.

- The Push Aequi flex is worn over the sock
- Place the brace around the foot. The rigid part is on the inside of the foot. Put the long straps to the sides
- Fasten the diagonal strap marked with  (fig. 1)
- Take hold of the elastic strap marked with  at the outside of the foot. Pull this diagonally across the instep, wrap it around the ankle just above the heel and fasten in place (fig. 2)
- Finally take hold of the elastic strap on the inside of the foot, marked with . Pull this from the instep across the dorsum (top of the foot), wrap it around the ankle over the first strap and fasten in place (fig. 3). The instep is the lower part of the foot, where the strap starts.

WASHING INSTRUCTIONS

The brace can be washed by hand at temperatures up to 30°C. Do not use any bleaching agents.

Hang the brace up to dry in the open air (not near a heating source or in a tumble drier).

To extend the brace's service life, it should be dried after use, if necessary, before being stored.

IMPORTANT

The efficient use of Push Aequi flex depends on the correct diagnosis of the reasons for pain or instability. Always check first with your doctor, a physiotherapist or another specialist before using a brace.

If this brace is to function effectively, it is important:

- to follow the washing instructions
- that the elastic material can still be stretched
- that the Velcro fastens efficiently
- that the brace is not severely damaged (no tears, loose seams)
- that the brace is not subjected to temperatures higher than 70°C.
- to please be aware that an ill-fitting brace does not function properly.

Push® is a registered trademark of Nea International bv
P.O. Box 305, 6199 ZN Maastricht-Airport, Netherlands
www.push.eu

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